


Experiment Report

In-situ PL Data



Time range: 0.2 - 37.6 s
Wavelength range: 40 - 60 nm

AI Agent Analysis

Analysis of Time-Resolved Photoluminescence Data

The experiment collected time-resolved photoluminescence (PL) spectra over a wavelength range of 40-60 nm and a time range of 0.2-37.6 seconds. The data shows several key features:

1. **Temporal Evolution:** The PL intensity generally decreases with time, indicating carrier recombination. The decay appears non-exponential, suggesting multiple recombination pathways or trap states.
2. **Spectral Features:** The PL spectrum shows a broad emission band with peak intensity around 50-52 nm. The spectral shape remains relatively constant over time, indicating stable emission characteristics without significant spectral diffusion.
3. **Intensity Distribution:** The highest intensities are observed at early time points (0.2-2 seconds) and decrease by approximately one order of magnitude over the measured time window. This suggests relatively long-lived excited states or slow recombination kinetics.
4. **Data Quality:** The signal-to-noise ratio is good at early times but decreases at longer times due to lower intensity counts. The repeated wavelength measurements (5 repeats per wavelength) provide good statistical confidence.

Recommended Next Steps:

1. Perform multi-exponential fitting of intensity decay at the peak wavelength to extract lifetime components.
2. Analyze spectral moments (center wavelength, full width at half maximum) as functions of time to detect spectral shifts or broadening.
3. Consider measuring at longer time scales (up to several minutes) to capture complete decay kinetics.
4. Repeat experiment under different excitation conditions (power, wavelength) to study excitation-dependent effects.
5. Perform temperature-dependent measurements to investigate thermal activation of non-radiative pathways.

The dataset provides a solid foundation for understanding the carrier dynamics in the studied material system. Further analysis with appropriate fitting models will yield quantitative parameters such as recombination lifetimes and quantum efficiency estimates.